



The IA100 uses a variety of industrial valves, flow meters and sensors, together with an OMRON PLC to demonstrate the principles of PID control.

The student is able to change the following parameters:

Set Value,

Proportional bandwidth,

Integral time,

Derivative time,

Sampling period and

Input filter coefficient

and observe the effects upon the system.

Manufacturer : **AMRON**

Tel : (033) 3868309

Fax : (033) 3868556

E MAIL: amron@pixie.co.za

P. O. BOX 2435
Pietermaritzburg
3201